INSTRUCTIONS FOR USE

Cavex Glass Ionomer Cement

A glass ionomer cement for cementation. Its outstanding features include low solubility in the mouth and low acidity.

Fields of application

For cementation of crowns, bridges, inlays, onlays, pins, posts and orthodontic bands.

Application

Clean enamel, dentine and metal surfaces thoroughly and dry completely. Make sure not to leave any traces of other materials, which could impair the properties of the glass ionomer luting cement. Temporary work should be attached preferably with an eugenol-free cement. Before attachment clean the core, rinse (e.g. with H₂O₂) and dry.

Shake powder well before use!

The working temperature should be 15–23°C.

Use one drop of liquid for one level measuring scoop of powder (powder / liquid ratio = 3.0:1 g/g)

Put powder and liquid on a suitable glass plate or a suitable mixing pad. When dispensing the liquid, hold the dropper bottle vertically.

It is recommended to divide the powder in 3 portions and to mix it portion by portion into the liquid, using a solid plastic spatula or another abrasion-resistant spatula.

The mixing time is approx. 30 seconds. The glass ionomer luting cement can be worked for approx. 3 minutes at room temperature and sets within 5–7 minutes (within 3–5 minutes in the mouth). Higher temperatures shorten the working time and lower temperatures lengthen it.

The correct consistency is achieved when the glass ionomer luting cement is applicable to the core with a hard brush. It is essential to avoid contact with water and saliva during application and setting.

Indications / Precautions

Excess material can be removed after 10–15 minutes, preferably with an explorer. Desensitize strongly ground and sensitive cores prior to a temporary attachment with a fluoride varnish. Remaining varnish and other substances (e.g. eugenol, isolating agent) can affect the adhesion and should be removed completely before using the luting cement.

Do not apply the cement excessively and maintain sufficient dentine thickness to avoid pulp irritations. If necessary, pulp protection (e.g. calcium hydroxide) should be applied in the proximity of the pulp or a core build-up should be carried out, if indicated.

Rinse instruments with cold water immediately after use and before the setting of the cement residues. After removal of orthodontic bands, cement residues can be removed gently and rapidly with an ultrasonic scaler.

The glass ionomer luting cement contains polyacrylic acid, fluoro-silicate and parabens. Do not use in cases of hypersensitivities against these ingredients.

Storage information

Close powder container immediately after use, in order to avoid infiltration of humidity, which may result in a loss of compressive strength.

Store at temperatures 4–23°C. Do not use after expiry date.

Our technical advice, whether verbal, in writing or by way of trials, is given in good faith but without warranty, and this also applies where proprietary rights of third parties are involved. It does not release you from the obligation to test the products supplied by us as to their suitability for the intended processes and uses. The application, use and processing of the products are beyond our control and, therefore, entirely your own responsibility. Should, in spite of this liability be established for any damage, it will be limited to the value of the goods delivered by us and used by you. We will, of course, provide consistent quality of our products within the scope of our General Conditions of Sale and Delivery.

CAVEX HOLLAND BV
Fustweg 5, 2031 CJ Haarlem, The Netherlands
T+31 23 530 77 00, dental@cavex.nl, www.cavex.nl

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