**Cavex SiliconA Bite**

Scannable polyvinyl siloxane for bite registration, superior hardness, very short time in mouth.

Cavex SiliconA Bite is an addition cross-linking polyvinyl siloxane with a very high final Shore D hardness of 40. Its surface properties and optimized color enable it to be optically recorded in CAD/CAM systems to portray the opposing dentition.

**Technical data**

Color	Working time** [min:s]	Minimum time in mouth** [min:s]	Shore D hardness	Shore A hardness
Lavender	0:30	0:30	40	≥ 94

\*1 Total working time (including mixing) at 23°C/73°F; shortened to 20s at oral temperature (35°C/95°F)  
\*2 Minimum time that the bite registration material must be in a closed mouth.

The data given are furnished to the best of our knowledge and are based on internal test results.

**Praktiske tips**  
Please visually inspect for damage before use. Do not use damaged products. Mixing tips and bite registration tips are for single use only. Dispense a small amount of material from a 50 ml cartridge before attaching a mixing tip. After use, leave the mixing tip on the cartridge to keep it closed. Immediately before the next use, remove the mixing tip and attach a new one.  
We use registration materials based on polyvinyl siloxane as touchless gloves, retraction solutions or freshly inserted composite fillings. It can change the curing behavior. Cavex SiliconA Bite should be at room temperature when used (23°C/73°F). Lower temperatures, for example when the materials are stored in a refrigerator can extend the working time, and higher temperatures can shorten it.  
Hardened polyvinyl siloxanes are chemically resistant. Clothes can be permanently stained if the contact material is not completely cured. Cavex SiliconA Bite is used with the Cavex SiliconA Mixing Tip green (1:1). It can also be combined with a bite registration tip.

**Uses**

- Standard bite registration in the end bite position
- Portrayal of opposing dentition
- Coating of bite forks
- Key material for needle point registration
- Production of small model segments

**Warning and safety instructions**

If the material contacts your eyes, rinse with copious amounts of water for several minutes. If symptoms continue, consult a doctor. Do not inhale, swallow or ingest. If health problems arise after swallowing impression material, immediately contact a doctor. In rare cases, intestinal blockage may arise. Immediately consult a doctor if the product is inhaled.  
Do not leave any residual material in the sulcus or oral cavity.  
Cavex SiliconA Bite should not be used for conventional impressions since its high final hardness makes it very difficult to remove from the mouth. Undesired tooth extractions may result.  
Hypersensitive persons may experience an allergic reaction to the product.  
Only to be used by qualified dental personnel.

**Disinfection with Cavex ImpreSafe**

Cavex SiliconA Bite can be disinfected by using Cavex ImpreSafe. Cavex ImpreSafe is a non-toxic, fast, safe, and very effective disinfectant. Three minutes are needed to kill bacteria, fungi and viruses. Due to this short contact time the surface of the impression remains 100% intact and negative impact on the end result is avoided.  
Cavex SiliconA Bite can also be treated with an aqueous sodium hypochlorite solution (5,25% NaOCl) for up to 10 minutes. The treatment does not affect the surface or dimensions. When using other disinfectant solutions suitable for silicones, follow the disinfection and processing instructions of the manufacturer.

**Storage conditions**

Do not store above 25°C/77°F

Do not use after the expiration date (see the information on the package).

Our technical advice, whether verbal, in writing or by way of trials, is given in good faith but without warranty, and this also applies where proprietary rights of third parties are involved. It does not release you from the obligation to test the products supplied by us as to their suitability for the intended processes and uses.  
The application, use and processing of the products are beyond our control and, therefore, entirely your own responsibility. Should, in spite of this liability be established for any damage, it will be limited to the value of the goods supplied by us and used by you. This also concerns a contractual, contractual or other legal relationship of our company with you under conditions of Sale and Delivery.

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