INSTRUCTIONS FOR USE

Cavex Non Gamma-2

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The amalgam made of Cavex Non Gamma-2 is free from the corrosion-prone and weak gamma-2 phase. Durable restorations with a lasting lustre can be made.

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Indication
Stress-bearing fillings in posterior teeth (Class I and II).

Contraindications
L'usage d'amalgam n'est pas recommandé when a physician has deemed it unsuitable for the patient.

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• Hypersensitivity (allergic reaction) to amalgam components. In such cases the use of suitable alternative materials must be considered.

Warnings
After placement or removal of dental amalgam restorations increased mercury concentration in blood and urine has been observed. According to available scientific knowledge this increase has not been associated with any adverse health effects.

• If a patient will be exposed to dental restorations containing mercury vapour or dental amalgam particles during handling, placement or removal of dental amalgam should be avoided. Placement and removal of dental amalgam fillings should be performed with appropriate water spray and vacuum suction.

• Instructions and regulations for storage and disposal of dental amalgam waste must be observed.

Directions for use

Cavex Non Gamma-2 is available as pre-dosed capsules. It should be noted, that the working and hardening time would be somewhat different, depending on the type of mixer used: a high-energy mixer gives a faster hardening. The mixed amalgam should be a plastic, coherent mass with a shiny surface.

Mixing ratio for Cavex Non Gamma-2 is: 10 parts alloy / 11 parts mercury by weight

Cavex Non Gamma-2 in all cases 5 - 7 seconds in High-energy mixer (such as Silamat)

- I-spill (alloy/mercury 400/435 mg)

- II-spill (alloy/mercury 600/653 mg)

- III-spill (alloy/mercury 800/865 mg)

Notes:

- 1. For safety reasons, it is best to replace the empty capsule after taking out the mixed amalgam, before cleaning or disposal.
- 2. The pre-dosed capsules (Type I) are self-curing; contact between alloy and mercury occurs automatically when the capsule is broken due to the action of the mercury. Length 32 mm, max. outer diameter: 14 mm.
- 3. Storage at temperatures no higher than 25°C is recommended.
- 4. Dispose of cotton wool and used capsules in accordance with national regulations.

Insertion and condensation
After mixing, the amalgam should immediately be inserted into the cavity using an amalgam carrier. The filling should be built up in layers under adequate condensation pressure using a plugger appropriate for the size of the cavity. Do not use an ultrasonic condenser. Carefully remove excess mercury.

Finishing and polishing
Approx. 8-10 minutes after mixing, the filling can be finished. Polishing should be done in a subsequent session.

"THIS ALLOY CONTAINS ZINC: THE AMALGAM MADE THEREFROM MAY SHOW EXCESSIVE EXPANSION IF MOISTURE IS INTRODUCED DURING MIXING OR COMPACTING"

Labelling
This information is in compliance with ISO24234 and also with European legislation.

MERCURY DANGER

- May impair fertility or cause harm to an unborn child.
- Lethal if inhaled.
- Causes damage to organs in case of long-term or repeated exposure.
- Highly toxic to aquatic organisms, may cause long-term adverse effects.

Recommendations for safety precautions
Dental clinical personnel, when using silver amalgam, should exercise proper procedures to avoid contact with the alloy. The use of a protective cap to cover the capsule will prevent the release of mercury vapour into the dental environment. Waste amalgam should be collected for disposal in a manner, which will protect both those who handle the waste and the environment.

1. Treat mercury and amalgam with care. Be MERCURY-AWARE.
2. Ideally a high-energy mixer should be equipped with a protective cap to cover the capsule while mixing.
3. Spills of mercury can be collected with the aid of a "mercury-collecting forceps" or a small amount of freshly mixed amalgam that will easily absorb liquid mercury. NEVER use a vacuum cleaner.

4. Removal of amalgam and fillings must be done under water-cooling and with effective vacuum suction. The use of eye protection and mouth mask is advisable.

5. The wastewater drain of the spittoon should be equipped with an amalgam separator. **Le déversement des déchets de l'amalgam waste should be taken into consideration.**

6. Excess mixed amalgam should be kept in a container, designed for the purpose, or at least in a closed container filled with photographic fixer.

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