QUADRANT TOTAL ETCH

(DE)	Gebrauchsanweisung	3
(GB)	Instructions for use	7
(FR)	Mode d'emploi	11
(ES)	Instrucciones de uso	15
(IT)	Istruzioni per l'uso	19
(NL)	Gebruiksaanwijzing	23
(GR)	Οδηγίες χρήσης	27
(CZ)	Návod k použití	31
(HU)	Használati utasítás	35
(LT)	Naudojimo instrukcija	39
(HR)	Upute za uporabu	43
(RO)	Instrucțiuni de utilizare	47
(TR)	Kullanma talimatı	51
(AE)	تعليمات الاستخدام	55
(BG)	Инструкции за употреба	59
(EE)	Kasutusjuhend	63
(SI)	Navodila za uporabo	67

CAVEX

Pastös
Pasty
Pâteux
Pastoso
Pastoso
Pasteus
πηλώδης
Pastovité
Pasztaszerű anyag
Pastos pavidalo
Pastozan
Pastă
Macun gibi
مادة عجينية
пастьозно
Pastakujuline
Pastozen

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QUADRANT TOTAL ETCH

Instructions for use 

Medical Device – for use by dental health care professionals only.

Intended purpose

Quadrant Total Etch is an etching agent based on orthophosphoric acid. It is used to etch enamel and condition dentine in the context of adhesive filling therapy.

Indications

Structural defects, alteration of tooth structure or missing teeth

Dental applications:

Enamel etching and enamel/dentine etching in the context of the total etch technique prior to the adhesive bonding of:

- Composite and compomer restorations
- Laboratory-fabricated, indirect restorations (e.g. inlays, crowns, bridges, veneers)
- Sealants

Contraindications

Application to dentine close to the pulp (less than 1 mm). In case of known or suspected allergy to components of the product, its use is contraindicated.

Patient target group

This dental material is suitable for patients requiring dental treatment for the mentioned indications with consideration of the contraindications. For pregnant and nursing women, the treatment risks must be weighed carefully against the benefits taking into consideration the unborn child or infant.

Intended user

For use by dental health care professionals only.

Clinical benefit

Cavex medical devices ensure the rehabilitation of oral functions such as chewing, speaking and aesthetics. They stabilize the remaining dentition and/or the alveolar ridge restoratively or prosthetically.

Potential side effects

Hypersensitivities to the product or its components cannot be excluded in individual cases.

Composition

Quadrant Total Etch contains: water, phosphoric acid, silicon dioxide.

Product application

Prior to treatment, carefully clean tooth with fluoride-free polishing paste. The use of a rubber dam is recommended. Visually inspect for damage before use. Damaged products must not be used.

1. After preparation in accordance with adhesive filling therapy, rinse the cavity and dry with oil-free air.
2. Starting from the bevelled enamel areas, apply Quadrant Total Etch to the entire cavity surface, including the dentine. Ensure that the enamel is conditioned for at least 20 s, a maximum of 30 s and the dentine for 15 s. If desired, only apply Quadrant Total Etch to the enamel and let set for 20 – 30 s. Etch unprepared enamel, e.g. in case of fissure sealing, for at least 30 s and no more than 60 s. In case of fluoridated enamel, etching for no more than 60 s may also be required.

Caution: Areas near the pulp should be covered with a suitable pulp protection material according to the manufacturer's instructions.

3. Use a water sprayer to completely rinse off Quadrant Total Etch for 20 s under suction. The cavity must be completely free of any remnants.

- 
4. Dry the cavity with oil-free air (do not dry out dentine!). The etched enamel surface should appear chalky white. In case of contamination of the tooth surface (e.g. with saliva), the etching must be repeated (for a maximum of 10 s dwell time).
 5. Perform further treatment with adhesive (e.g. Quadrant Bondings) and composite (e.g. Quadrant composites) in accordance with manufacturer directions.

Quadrant Total Etch may be applied using a brush or cannula. Both brushes and cannulas are for one-time use only.

Release the plunger after applying the gel. When using syringes in mouth, use common hygienic protective sleeves to prevent contamination. Do not reuse after contamination. After use, close Quadrant Total Etch with the appropriate plug or cap.

Interactions

Liner and base materials, such as calcium hydroxide, may be dissolved by the acid, depending on their type and composition.

General information

Please pay attention to the following safety instructions and all processing instructions in other sections. This product must be used according to the instruction for use only. Any use not complying with these instructions is at the discretion and sole responsibility of the dental professional itself.

Warnings and safety instructions

Quadrant Total Etch contains 35 % phosphoric acid and can cause chemical burns. May cause skin and eye irritation. Avoid contact with skin. May cause an allergic skin reaction. Wear suitable protective gloves. In case of contact with skin, immediately wash with plenty of water and soap. Avoid contact with eyes. Eye

protection is recommended for the patient. In case of contact with eyes, rinse with water for several minutes. If symptoms persist, consult an eye specialist.

The etching agent must not be swallowed by the patient, otherwise incompatibilities may occur.

Special advice

According to EU Medical Device Regulation users/ patients are obliged to report serious incidents with a medical device to the manufacturer and to the competent authority of the country, the incident occurred.

Storage conditions

Store below 25 °C (77 °F). Protect from sunlight. Store in a dry place. Do not use after expiry date.

Disposal information

Recommendation: Dispose of in accordance with official regulations. Do not dispose of contents or partially emptied packaging in the household waste or allow them to enter the sewage system.

European Waste Catalogue: 180106 Chemicals consisting of or containing hazardous substances.

Complaint notice

Please quote batch number and article number or UDI in all correspondence about the product.

Article number **REF** and batch number **LOT** or UDI: See note on packaging.

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Kulzer GmbH