QUADRANT UNI-1-BOND

(DE) Gebrauchsanweisung	3	(LT) Naudojimo instrukcija	68
(GB) Instructions for use	10	(HR) Upute za uporabu	75
(FR) Mode d'emploi	17	(RO) Instrucțiuni de utilizare	82
(ES) Instrucciones de uso	24	(TR) Kullanma talimatı	89
(IT) Istruzioni per l'uso	32	(AE) تعليمات الاستخدام	96
(NL) Gebruiksaanwijzing	39	(BG) Инструкции за употреба	102
(GR) Οδηγίες χρήσης	46	(EE) Kasutusjuhend	110
(CZ) Návod k použití	54	(SI) Navodila za uporabo	117
(HU) Használati utasítás	61		

66112771/01

CAVEX

L2605_27633_A10012_66112771_GB Quadrant UNI-1-BOND e-IFU.indd 1

05.05.26 10:4

KULZER

MITSUBI CHEMICALS GROUP

Dateiname: KUL2605_27633_A10012_66112771_GB Quadrant Uni 1 Bond e-IFU		ActiVote-Nr: A10012	Maße: 210 x 74 mm	Faltmaß: 105 x 74 mm (geheftet)
SAP-Nr: 66112771	Version: 01	AK-Version: Print-PDF	Article: PB Quadrant Uni 1 Bond e-IFU	Datum: 05.05.2026
Projektmanager: Markus Weik		Druckfarben: Schwarz	Freigabe am: Offset	

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Stand: 2026-04

QUADRANT UNI-1-BOND Instructions for use

Medical Device – for use by dental health care professionals only.

Intended purpose

Quadrant Uni-1-Bond is a light-cured single-component adhesive for use in adhesive restorative dentistry. Quadrant Uni-1-Bond has been developed for the adhesive bonding of plastic filling materials (e.g. composites, compomer) to dental hard substances and laboratory-fabricated restorations (e.g. ceramics). With Quadrant Uni-1-Bond, priming, bonding and desensitizing can be done in one step. Quadrant Uni-1-Bond is a light-activated adhesive in an ethanolic solution. Before using Quadrant Uni-1-Bond, the dental hard substance is conditioned with an etching gel (etch&rinse).

Indications

Structural defects, alteration of tooth structure or missing teeth

Dental applications:

- Adhesive fixing of direct composite and compomer restorations
- Adhesive fixing of indirect laboratory-fabricated ceramic and composite restorations (inlays, onlays, veneers and crowns)
- Treatment of hypersensitive tooth regions

Contraindications

In case of known or suspected allergy to components of the product, its use is contraindicated. Do not use this product in case of known or suspected allergies against (meth)acrylate compounds or glutaraldehyde. Do not bring Quadrant Uni-1-Bond into direct contact with opened pulp tissue.

Patient target group

This dental material is suitable for patients requiring dental treatment for the mentioned indications with consideration of the contraindications. For pregnant and nursing women, the treatment risks must be weighed carefully against the benefits taking into consideration the unborn child or infant.

Intended user

For use by dental health care professionals only.

Clinical benefit

Cavex medical devices ensure the rehabilitation of oral functions such as chewing, speaking and aesthetics. They stabilize the remaining dentition and/or the alveolar ridge restoratively or prosthetically.

Composition

Quadrant Uni-1-Bond contains: UDMA, ethanol, HEMA, water, methacrylated polycarboxylic acid, 4-META, glutaraldehyde

Potential side effects

Hypersensitivities to the product or its components cannot be excluded in individual cases.

Product features

Available in a 4 ml bottle

Application

1. Adhesive fixing of direct composite and compomer restorations

1.1 Preparation

- Clean tooth with a toothpaste free from oil and fluoride.
- Prepare cavity as specified by the rules of adhesive filling treatment.
- Rinse prepared site with water and blow-dry. Avoid surface contamination.

- Areas near the pulp should be covered with a suitable pulp protection material according to the manufacturer's instructions. Cover the dentine near the pulp only. Leave the remainder of the cavity free for adhesion of Quadrant Uni-1-Bond.

1.2 Conditioning of the tooth surface

Before applying Quadrant Uni-1-Bond, the prepared cavity must be pretreated with a suitable etching gel according to the manufacturer's instructions (e.g. Quadrant Total Etch from Cavex) and thoroughly rinsed with water and dried. Overdrying of the dentin should be avoided; the dentin surface should be slightly damp.

1.3 Application of Quadrant Uni-1-Bond

- Dispense Quadrant Uni-1-Bond in a mixing well, extract it with an applicator tip or a soft disposable brush and apply to the entire surface of the cavity.
- Allow Quadrant Uni-1-Bond to soak in for 15 seconds.
- Quadrant Uni-1-Bond contains ethanol as a solvent. Carefully evaporate solvent and residual moisture under a gentle air stream until no movement of liquid can be detected.

An excessively strong air flow at the start of blow-drying will dilute the bonding and may result in inadequate adhesion.

- The surface must be visibly shiny after application of Quadrant Uni-1-Bond and after evaporation of the solvent. Make certain that the complete surface of the cavity is fully covered. If the cavity does not appear to be completely shiny all over, apply Quadrant Uni-1-Bond a second time as described above. Make sure that no excess remains in corners of the cavity. Cure Quadrant Uni-1-Bond for 20 seconds using a suitable halogen or LED polymerisation lamp (wavelength 440–470 nm, light intensity at least 600 mW/cm²). If the light power is too low there will be insufficient adhesion.

Insufficient lamp intensity will result in inadequate adhesion. Polymerisation lamps should be tested with reliable testing devices at regular intervals. The light emitting window should be placed as close as possible to the plastic surface during polymerisation.

- Then immediately apply and process the filling material in the cavity as directed by the manufacturer.

2. Adhesive fixing of indirect laboratory-fabricated ceramic and composite restorations (inlays, onlays, veneers and crowns)

2.1 Preparation of the restoration

- Prepare restoration surface as directed by the manufacturer's instructions for use.

2.2 Conditioning of the tooth surface

Proceed as described in 1.2

2.3 Application of Quadrant Uni-1-Bond

Make sure that no excess material is left in the corners of the cavity and that the Quadrant Uni-1-Bond coating is carefully blow-dried before polymerisation.

Otherwise proceed as described in 1.

- Polymerise Quadrant Uni-1-Bond for 20 s with a halogen or LED polymerisation lamp (see instructions above).
- Place composite cement (self-curing or dual-curing) according to the manufacturer's instructions for use.

3. Treatment of hypersensitive tooth regions and cleaning the tooth

- Carefully rinse the surface with water and blow-dry.

3.1 Conditioning hypersensitive tooth regions

Proceed as described in 1.

3.2 Application of Quadrant Uni-1-Bond

Proceed as described in 1.

- Polymerise Quadrant Uni-1-Bond for 20 s with a halogen or LED polymerisation lamp (see instructions above).
- Carefully remove the oxygen inhibition layer with an alcohol-soaked pellet.

If Quadrant Uni-1-Bond does not sufficiently desensitise, repeat the process as described above.

General information

Please pay attention to the following safety instructions and all processing instructions in other sections. This product must be used according to the instruction for use only. Any use not complying with these instructions is at the discretion and sole responsibility of the dental professional itself.

- Visually inspect for damage before use. Damaged products must not be used.
- Materials containing eugenol may impair the polymerisation of Quadrant Uni-1-Bond.
- Ensure product is at room temperature prior to application.

Warnings and safety instructions

Highly flammable liquid and vapour. May cause skin and eye irritation. Avoid contact with skin. May cause an allergic skin reaction. Wear suitable protective gloves. In case of contact with skin, immediately wash with plenty of water and soap. Avoid contact with eyes. In case of contact with eyes, rinse with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If symptoms persist, consult an eye specialist. Toxic to aquatic life with long lasting effects.

- Quadrant Uni-1-Bond contains ethanol. Ethanol is very flammable. Do not inhale vapour.

- Quadrant Uni-1-Bond contains glutaraldehyde. If it is inhaled or swallowed, it may damage health. Seek medical attention immediately if swallowed. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources - no smoking. May cause respiratory irritation. Avoid breathing dust/fume/gas/mist/vapours/spray. Avoid static charging. Keep away from all sources of ignition.
- Use only in sufficiently well-ventilated rooms. Do not inhale vapours.
- Do not apply adhesive directly to the exposed pulp.
- Eye protection is recommended for the patient. The use of a rubber dam is strongly recommended.
- Close the Quadrant Uni-1-Bond bottle carefully after each use.
- To prevent cross-contamination, do not reuse bottles after contamination with bodily fluids.
- Saliva or blood in the field of application during treatment may lead to failure. Use rubber dam or other suitable means of isolation.

Special advice

For SSCP please see EUDAMED database (<https://ec.europa.eu/tools/eudamed>) or contact sscp-service@kulzer-dental.com

To find the SSCP in the EUDAMED use the UDI code attached at the packaging of the product.

According to EU Medical Device Regulation users/ patients are obliged to report serious incidents with a medical device to the manufacturer and to the competent authority of the country, the incident occurred.

Storage conditions

Do not use the material after the expiration date. After opening, use by expiry date. Store at temperatures of -10 to 25°C (14 to 77°F). The product can be stored in refrigerator. Use the adhesive at room temperature only. Do not expose to direct sunlight. Store in a sufficiently well ventilated dry place. The shelf life will be reduced by unsuitable storage conditions and the function of the product may be affected. Close the container tightly immediately after use.

Disposal information

Recommendation: Dispose of in accordance with official regulations. Do not dispose of contents or partially emptied packaging in the household waste or allow them to enter the sewage system.

European Waste Catalogue: 180106 Chemicals consisting of or containing hazardous substances.

Complaint notice

Please quote batch number and article number or UDI in all correspondence about the product.

Article number **REF** and batch number **LOT** or UDI: See note on bottle/packaging.

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Kulzer GmbH