QUADRANT UNI-SE-BOND

(DE) Gebrauchsanweisung	3	(LT) Naudojimo instrukcija	57
(GB) Instructions for use	9	(HR) Upute za uporabu	63
(FR) Mode d'emploi	15	(RO) Instrucțiuni de utilizare	69
(ES) Instrucciones de uso	21	(TR) Kullanma talimatı	75
(IT) Istruzioni per l'uso	27	(AE) تعليمات الاستخدام	81
(NL) Gebruiksaanwijzing	33	(BG) Инструкции за употреба	87
(GR) Οδηγίες χρήσης	39	(EE) Kasutusjuhend	93
(CZ) Návod k použití	45	(SI) Navodila za uporabo	99
(HU) Használati utasítás	51		

66112772/01

CAVEX

L2605_27623_A10010_66112772_GB Quadrant Uni-Se-Bond e-IFU.indd 1

06.05.26 11:00

KULZER
wpsuj orthopaedics group

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QUADRANT UNI-SE-BOND

Instructions for use (GB)

Medical Device – for use by dental health care professionals only.

Intended purpose

Quadrant Uni-SE-Bond is a light-curing, self-conditioning all-in-one adhesive for use in adhesive, restorative dentistry. Quadrant Uni-SE-Bond was developed for the adhesive bonding of direct filling materials.

As a universal adhesive, it may be used as follows depending on the indication:

- SELF-ETCH TECHNIQUE for etching, priming, bonding and desensitisation in a single step.
- SELECTIVE ETCH TECHNIQUE: adhesion to the tooth enamel can be maximised by applying selective enamel etching with phosphoric acid.
- ETCH & RINSE TECHNIQUE (total-etch) with an additional etching step with phosphoric acid on enamel and dentine.

Indications

Structural defects, alteration of tooth structure or missing teeth

Dental applications:

1. Bonding of direct restorations for all cavity classes (Black) using light-curing, dual-curing or self-curing methacrylate based composites/compomers.
2. Bonding of light-curing, dual-curing or self-curing core build-up materials.
3. Sealing of hypersensitive tooth areas.
4. Sealing of cavities prior to amalgam restorations.
5. Bonding of fissure sealants.
6. Sealing of cavities and core preparations prior to temporary cementation of indirect restorations (according to the immediate dentine sealing technique).

Contraindications

In case of known or suspected allergy to components of the product, its use is contraindicated. Do not use this product in case of known or suspected allergies against (meth)acrylate compounds. Do not allow Quadrant Uni-SE-Bond to come into direct contact with exposed pulp tissue.

Patient target group

This dental material is suitable for patients requiring dental treatment for the mentioned indications with consideration of the contraindications. For pregnant and nursing women, the treatment risks must be weighed carefully against the benefits taking into consideration the unborn child or infant.

Intended user

For use by dental health care professionals only.

Clinical benefit

Cavex medical devices ensure the rehabilitation of oral functions such as chewing, speaking and aesthetics. They stabilize the remaining dentition and/or the alveolar ridge restoratively or prosthetically.

Potential side effects

This product or one of its components may, in particular cases, cause hypersensitive reactions. In the case of a suspected allergy, we recommend taking an allergy test prior to use.

Composition

Acetone, UDMA, water, 4-META, MDP

Application steps

Preparation of the tooth surface (except fissure sealants)

- Prepare cavity according to the rules of the adhesive technique.
- Clean with water spray and dry with water-free and oil-free air.
- Areas near the pulp should be covered with a suitable pulp protection material according to the manufacturer's instructions.

Application options

Quadrant Uni-SE-Bond may be used as a self-etching adhesive in the self-etch technique.

For selective enamel etching, apply 35 % phosphoric acid etching gel (e.g. Quadrant Total Etch) to the prepared enamel and allow to work for 20–30 seconds. Then rinse thoroughly with water and dry with an oil-free air flow. Do not overdry.

In the etch & rinse technique, apply 35 % phosphoric acid etching gel (e.g. Quadrant Total Etch) to the prepared enamel and dentine and allow to work (20–30 seconds on enamel/15 seconds on dentine). Then rinse thoroughly with water and dry with an oil-free air flow. Do not overdry.

- Fissure sealing: after pretreating with polishing paste or powder jets and water spray, etch with 35 % phosphoric acid etching gel (e.g. Quadrant Total Etch) for 20–30 seconds, rinse thoroughly, dry and apply adhesive.

Caution: Pay attention to the instructions for using the etching gel.

Preparation of Quadrant Uni-SE-Bond

From the bottle:

- Drop Quadrant Uni-SE-Bond into the mixing well and use immediately (within 3 minutes). When product has been removed, reseal carefully.

Application to tooth surface

- Use the applicator tip or brush to apply Quadrant Uni-SE-Bond to the entire surface of the cavity.
- Gently rub in the adhesive for 20 seconds.
- Air dry carefully with a gentle oil-free air flow until the adhesive film no longer moves. Move the air flow from outside to inside with increasing intensity at the same time. Caution: If the air flow at the start of blowing is too strong, this leads to the bonding agent thinning out and may result in poor adhesion and post-operative discomfort. Do not leave any excess bonding material/puddles, e.g. on the cavity base.
- The surface must be visibly glossy both after applying. If the cavity surface does not appear universally glossy, apply Quadrant Uni-SE-Bond again as described. Avoid surface contamination.
- Polymerise Quadrant Uni-SE-Bond for 10 seconds with a conventional curing lamp (wavelength range 440–480 nm; light intensity min. 600 mW/cm²).
- Introduce restoration material immediately in accordance with the manufacturer's instructions.
- Caution: If the light power is too low this leads to insufficient adhesion. Curing lamps should be inspected at regular intervals with calibrated test equipment. The light emission window should be placed as close as possible to the cavity surface during polymerisation.
- Caution: Carefully remove oxygen inhibited layer with an alcohol soaked pellet after light curing.

General information

Please pay attention to the following safety instructions and all processing instructions in other sections. This product must be used according to the instruction for use only. Any use not complying with these instructions is at the discretion and sole responsibility of the dental professional itself.

- Visually inspect for damage before use. Damaged products must not be used.
- Materials containing eugenol may impair the polymerisation of Quadrant Uni-SE-Bond.
- Quadrant Uni-SE-Bond is light curing and should be protected against ambient light.
- The use of disinfectants, astringents and rinsing solutions containing EDTA or hydrogen peroxide is not recommended. Residues of these materials may impair the adhesion and curing reaction of the adhesive. But Quadrant FiniSense and chlorhexidine solutions may be used.

Warnings and safety instructions

Highly flammable liquid and vapour. May cause skin and eye irritation. Avoid contact with skin. May cause an allergic skin reaction. Wear suitable protective gloves. In case of contact with skin, immediately wash with plenty of water and soap. Avoid contact with eyes. In case of contact with eyes, rinse with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If symptoms persist, consult an eye specialist. May cause respiratory irritation. Avoid breathing dust/fume/gas/mist/vapours/spray. May cause drowsiness or dizziness. Harmful to aquatic life with long lasting effects. The adhesive must not be swallowed by the patient, otherwise incompatibilities may arise.

- Eye protection is recommended for the patient. The use of a rubber dam is strongly recommended.
- Seal Quadrant Uni-SE-Bond bottles carefully after each use.
- Do not reuse bottles after contamination with body fluids in order to prevent cross-contamination.
- Saliva or blood in the field of application during treatment may lead to failure. Use rubber dam or other suitable means of isolation.

Special advice

For SSCP please see EUDAMED database (<https://ec.europa.eu/tools/eudamed>) or contact sscp-service@kulzer-dental.com
To find the SSCP in the EUDAMED use the UDI code attached at the packaging of the product.

Storage conditions

Do not expose to direct sunlight. Store in a sufficiently well ventilated dry place. Store at room temperature, not above 25°C (77°F). Incorrect storage may shorten the shelf-life and the product may prematurely lose its effectiveness. Do not use the material after the expiration date.

Disposal information

Recommendation: Dispose of in accordance with official regulations. Do not dispose of contents or partially emptied packaging in the household waste or allow them to enter the sewage system.
European Waste Catalogue: 180106 Chemicals consisting of or containing hazardous substances.

Complaint notice

Please quote batch number and article number or UDI in all correspondence about the product.
Article number ^(REF) and batch number ^(LOT) or UDI: See note on bottle/packaging.

According to EU Medical Device Regulation users/patients are obliged to report serious incidents with a medical device to the manufacturer and to the competent authority of the country, the incident occurred.

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- Als er zich tijdens de behandeling speeksel of bloed bevindt in het toepassingsgebied, kan dit leiden tot falen. Gebruik een rubberdam (bijv. ©) of andere geschikte isolatiemiddelen.
Kulzer GmbH