QUADRANT UNIVERSAL LC

GB	Gebruiksaanwijzing
GB	Instructions for use
FR	Mode d'emploi
ES	Instrucciones de uso
IT	Istruzioni per l'uso
NL	Gebruiksaanwijzing
GR	Οδηγίες χρήσης
CZ	Návod k použití
HU	Használati utasítás
LT	Naudojimo instrukcija

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QUADRANT UNIVERSAL LC Instructions for use

Medical Device – for use by dental health care professionals only.

Intended purpose
Quadrant Universal LC is a light curing, radio-opaque microhybrid composite that is using Microglas® technology. It is indicated for adhesively bonded, both colored restorations (anterior and posterior teeth) including occlusal surfaces. Its radiopacity is equivalent to 2.0–2.5 mm of dentin (anterior = 1 mm, enamel = 2 mm).

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Contraindications
Application to dentine should be no less than 1 mm. Use of this product is contraindicated in case of known or suspected allergic against components of this product. Do not use this product in case of known or suspected allergies against bisphenol-Acrylate-derivatives.

Patient target group
This dental material is suitable for patients requiring dental treatment for the mentioned indications with consideration of the mentioned indications. The use of this product is contraindicated in case of known or suspected allergies against mentioned components of this product. Do not use this product in case of known or suspected allergies against bisphenol-Acrylate-derivatives.

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Compositum

Barium Aluminium Boro Fluor Silicate Glas; BioGlas; Feldspat; TEGOMA, UDMA, Silicone Dioxide; Polymer; Titanium Dioxide; Fluorescent Pigments; Metalla Oxide Pigments; Organic Pigments; Aminobenzoatesolventes; BHT; Camphorone. Quadrant Universal LC contains approximately 61% filler by volume.

This material is intended for use by dental professionals. Quadrant Universal LC is intended for adhesive bonding of dental veneers, dental veneers with mirror-finish surfaces, including occlusal surfaces. Its radiopacity is equivalent to 2.0–2.5 mm of dentin (anterior = 1 mm, enamel = 2 mm).

Indications
Structural defects, alteration of tooth structure or missing teeth. Dental applications:

• Direct restorations of all cavity classes • Direct composite veneers • Shape corrections of teeth (i.e. diastemas, congenital defects in teeth, etc.) • Splicing of teeth • Restoration of primary teeth • Repairs of porcelain, composite (in combination with an adequate repair-system)

Contraindications
Application to dentine should be no less than 1 mm. Use of this product is contraindicated in case of known or suspected allergies against components of this product. Do not use this product in case of known or suspected allergies against bisphenol-A-derivatives.

Patient target group
This dental material is suitable for patients requiring dental treatment for the mentioned indications with consideration of the mentioned indications. The use of this product is contraindicated in case of known or suspected allergies against mentioned components of this product. Do not use this product in case of known or suspected allergies against bisphenol-A-derivatives.

Intended user
For use by dental health care professionals only.

Clinical benefit
Cave medical devices ensure the rehabilitation of oral functions such as chewing, speaking and aesthetics. They stabilize the remaining dentition and/or the avoidable ridge restorativity or prophylactically.

Potential side effects
This product or one of its components may in particular cause cause hypersensitive reactions. If suspected, information on the ingredients can be obtained from the manufacturer.

Indications
Structural defects, alteration of tooth structure or missing teeth. Dental applications:

• Direct restorations of all cavity classes • Direct composite veneers • Shape corrections of teeth (i.e. diastemas, congenital defects in teeth, etc.) • Splicing of teeth • Restoration of primary teeth • Repairs of porcelain, composite (in combination with an adequate repair-system)

Contraindications
Application to dentine should be no less than 1 mm. Use of this product is contraindicated in case of known or suspected allergies against components of this product. Do not use this product in case of known or suspected allergies against bisphenol-A-derivatives.

Patient target group
This dental material is suitable for patients requiring dental treatment for the mentioned indications with consideration of the mentioned indications. The use of this product is contraindicated in case of known or suspected allergies against mentioned components of this product. Do not use this product in case of known or suspected allergies against bisphenol-A-derivatives.

Intended user
For use by dental health care professionals only.

Clinical benefit
Cave medical devices ensure the rehabilitation of oral functions such as chewing, speaking and aesthetics. They stabilize the remaining dentition and/or the avoidable ridge restorativity or prophylactically.

Potential side effects
This product or one of its components may in particular cause cause hypersensitive reactions. If suspected, information on the ingredients can be obtained from the manufacturer.

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