SILICON A

CAVEX

(DE)	Gebrauchsanweisung	5
(GB)	Instructions for use	11
(FR)	Mode d'emploi	17
(ES)	Instrucciones de uso	23
(NL)	Gebruiksaanwijzing	29
(GR)	Οδηγίες χρήσης	35
(CZ)	Návod k použití	41
(HU)	Használati utasítás	47
(LT)	Naudojimo instrukcija	53
(HR)	Upute za uporabu	59
(RO)	Instrucțiuni de utilizare	65
(CN)	使用说明	71
(SK)	Návod na použitie	76
(BG)	Инструкции за употреба	82
(EE)	Udhëzime përdorimi	88
(RS)	Uputstva za upotrebu	94

KULZER

INTUBE CHEMICALS GROUP

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D = Mischzeit

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K = Intraoral Tip

Medical Device – for use by dental health care professionals only.

Intended purpose

Addition curing polyvinyl siloxane for impressions, for direct application.

Indications

Structural defects, alteration of tooth structure or missing teeth Dental applications: Indicated for all impression techniques.

Contraindications

The use of this product is contraindicated in case of known or suspected allergies against components of this product. The product contains allergens: peppermint oil.

Patient target group

This dental material is suitable for patients requiring dental treatment for the mentioned indications with consideration of the contraindications. For pregnant and nursing women, the treatment risks must be weighed carefully against the benefits taking into consideration the unborn child or infant.

Intended user

For use by dental health care professionals only.

Clinical benefit

Cavex medical devices ensure the rehabilitation of oral functions such as chewing, speaking and aesthetics. They stabilize the remaining dentition and/or the alveolar ridge restoratively or prosthetically or they aid in the rehabilitation process.

Potential side effects

Hypersensitivities to the product or its components cannot be excluded in individual cases. In case of suspicion please contact the manufacturer for ingredients information.

Composition

Elastomeric, silicone-based precision impression material: polyvinyl siloxane, addition curing

Material

SiliconA Putty
SiliconA Heavy Body
SiliconA Injection

Pack

Tub, 1:1, 2 x 250 ml
Cartridge, 1:1, 50 ml
Cartridge, 1:1, 50 ml

Product application

The SiliconA range is an optimally matched system and its different types of material can be combined for all impression techniques to suit specific requirements.

Not to be used in conjunction with other silicone materials. Note for implants: Follow the implant manufacturer's instructions regarding the impression technique and procedure steps.

Tubs must be carefully sealed after use. Do not inter-change caps/lids/spoons. Dispense and discard a small amount of material from a 50 ml cartridge before attaching a mixing tip. After use, leave the mixing tip on the cartridge to keep it closed. Immediately before the next use, remove the mixing tip, dispense and discard a small amount of material. Attach a new mixing tip, and again dispense and discard a small amount of material. Mixing tips and intraoral tips are single-use products. To prevent cross-contamination between patients, mixing tips and intraoral tips must be discarded after each use.

All products in the SiliconA assortment are dosed in a ratio of 1 : 1 ratio (volume or weight)

SiliconA	EN ISO 4823	Product colour		Mixing ratio	Mixing time	Extraoral working time* (23°C (73°F))	Intraoral working time (35°C (95°F))	Minimum time in mouth	Part of packaging material	
		Base	Catalyst						Mixing Tips	Intraoral Tips
SiliconA Putty	Type 0 putty	yellow	light grey	1:1	30	2:00	–	3:30	n.a.	n.a.
SiliconA Heavy Body	Type 2 medium bodied	green	light grey	1:1	auto	2:00	1:00	3:30	green AS95	white AS97
SiliconA Injection	Type 3 light bodied	blue	light grey	1:1	auto	2:00	1:00	3:30	yellow AS94	yellow AS96

*At room temperature of 23°C (73°F), 50% rel. humidity. Higher temperatures shorten, lower temperatures prolong these times.

In all addition-curing silicone impression materials, the setting process may be affected by contact with products containing active Sulfur or Nitrogen, or by fresh composite fillings or superstructures. These products include latex gloves, Polysulfide impression materials, retraction solutions, etc. SiliconA contains surfactants for hydrophilation, which in some patients may cause a superficial anesthetic effect.

Surplus material on the impression tray or in the patient's mouth should be avoided as this may result in gagging, swallowing or aspiration. To further reduce these risks during impression taking, the patient should be positioned as upright as possible.

Do not leave any residual material in the sulcus or oral cavity, as this can lead to abscess formation.

If there are substantial undercuts, highly divergent teeth, or other critical conditions in the impression area, appropriate blocking out must be performed when using SiliconA impression materials.

Model fabrication

The impression can be cast 60 minutes after removal from the mouth – there are no other time constraints. The impression is stable in long-term storage.

Recommended model material

All standard model materials can be used.

We recommend:

Class III plaster, MOLDANO®

Class IV plaster, MOLDASTONE®/MOLDASYNT®

Disinfection

Impressions can be treated with an aqueous sodium hypochlorite solution (5.25%, for up to 10 minutes*) or with other disinfectant solutions suitable for silicones. When doing so, please observe the disinfection/processing instructions issued by the relevant manufacturers.

** Gelson Luís Adabo, DDS, PhD, Elaine Zanarotti, Renata Garcia Fonseca, and Carlos Alberto dos Santos Cruz, DDS, PhD – Effect of disinfectant agents on dimensional stability of elastomeric impression materials. J Prosthet Dent, 81, 621-4, 1999*

Adhesion

We recommend using Kulzer Universal Adhesive.

Storage conditions

Do not store above 25°C (77°F). Store in a dry place. Do not use after expiry date.

Warnings and safety instructions

Polyvinyl siloxanes are chemically resistant. Material which has not set can soil clothing. **Wear protective gloves according to EN 420, category III, with protection against chemical risks according to EN 374, type C.** Avoid contact with skin and eyes. May cause skin and eye irritation. In case of contact with skin, immediately wash with plenty of water and soap. In case of contact with eyes, rinse with water for several minutes. If symptoms persist, consult an eye specialist. Do not swallow or ingest. If health problems arise after swallowing impression material, immediately contact a doctor. Intestinal blockage may arise in rare cases.

General information

Please pay attention to the following safety instructions and all processing instructions in other sections. This product must be used according to the instruction for use only. Any use not complying with these instructions is at the discretion and sole responsibility of the dental professional itself.

Application notes

Visually inspect for damage before use. Damaged products must not be used. SiliconA impression materials should be at room temperature 23°C (73°F) for processing.

Disposal information

Recommendation: Dispose of in accordance with official regulations. Do not dispose of contents or partially emptied packaging in the household waste or allow them to enter the sewage system.
European Waste Catalogue: 180106 Chemicals consisting of or containing hazardous substances.

Complaint notice

Please quote batch number and article number or UDI in all correspondence about the product. Article number **REF** and batch number **LOT** or UDI: See note on packaging.

According to EU Medical Device Regulation users/ patients are obliged to report serious incidents with a medical device to the manufacturer and to the competent authority of the country, the incident occurred.

Dated: 2026-05

Description of applicational / packaging icons

A = Only use pastes with same LOT No
C = Mixing ratio
D = Mixing time

E = Extraoral working time
F = Intraoral working time
G = Time in mouth
H = Static mixing tip
K = Intraoral tip

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Kulzer GmbH