



Cavex Holland BV
 P.O. Box 852
 2003 RW Haarlem
 The Netherlands

T +31 23 530 77 00
 F +31 23 535 64 82
 dental@cavex.nl
 www.cavex.nl

Cavex VacuFormer System Accessories

Cavex Instant Stone

Cavex Instant Stone is a unique plaster composite that casts perfect models in just 10 minutes.




article no. pieces
 VT010 72 pcs

Cavex Mouth Protector

Cavex mouth protector laminates are professional 3.8mm dual layer and 5mm triple layer laminates with a combined hard and soft layer.

Cavex Mouth Protectors (3,8mm)

colour	article no.	pieces
Red/White/Blue	VT020	12 pcs
Red	VT034	12 pcs
Blue	VT035	12 pcs
Clear	VT036	12 pcs



Cavex Mouth Protectors X-Pro (5mm)

colour	article no.	pieces
Black	VT032	12 pcs
Green	VT031	12 pcs
Clear	VT033	12 pcs

Cavex Splint

Cavex Splint laminates are made of a transparent, rigid material ideally suited for splints, surgical trays, orthodontic retainers, etc.

Cavex Splint (0,6mm)

colour	article no.	pieces
Transparent	VT021	50 pcs



Cavex Splint X-heavy (1,5mm)

colour	article no.	pieces
Transparent	VT026	25 pcs



Cavex Bleach

Cavex Bleach laminates are made of a soft/clear material that is easy to cut, resulting in a perfect finish. The soft material is very comfortable against the gum tissue and thus ideal for making bleaching spoons.

Cavex Bleach (1,0mm)

colour	article no.	pieces
Transparent	VT022	25 pcs



Cavex Bleach heavy (1,5mm)

colour	article no.	pieces
Transparent	VT030	25 pcs



Cavex Bruxism

Cavex Bruxism laminates are hard and crystal clear. They are 2mm thick. This plate is perfect for the treatment of bruxism and helps to prevent wear and tear of the elements.

Cavex Bruxism (2,0mm)

colour	article no.	pieces
Transparent	VT027	25 pcs



Cavex Bleach Tray Spacer

Cavex Bleach Tray Spacer is an easy-to-use blue coloured resin specially developed for creating reservoirs in bleaching spoons.

Cavex Bleach Tray Spacer

article no.	pieces
VT025	1 pcs



For more information, please visit www.cavex.nl or send an email to dental@cavex.nl

ZV912.02