

PREMIUM IMPRESSION MATERIALS

CAVEX SILICONA

-

A fully compatible, high quality system for every clinical situation.



YOUR
IMPRESSION
IS OUR
SPECIALTY

CAVEX

PREMIUM IMPRESSION MATERIALS

CAVEX SILICONA

Cavex has been a specialist in high-quality impression materials for over 70 years. Cavex SiliconA is an advanced silicone system with outstanding properties for all clinical situations.

Thanks to its superior hydrophilic properties, perfect impressions are possible, even in the presence of moisture. SiliconA is non-sticky, making it easy to handle. The impression can be poured 60 minutes after removal using any commercial plaster.

The material does not adhere to surfaces such as gloves, simplifying tray loading. Due to the high quality of the raw materials, impressions remain stable for at least 7 days and capture even the finest details.

Cavex SiliconA is available in three consistencies that complement each other perfectly:

- Yellow putty: kneadable, very high consistency
- Green Heavy Body (Monophase): medium consistency
- Blue Injection Type: low consistency

This makes both the sandwich and double impression techniques possible. Thanks to its thixotropic properties, the material only flows under pressure. The platinum catalyst system ensures optimal polymerization.

Features

- Hydrophilic and not sticky
- Long term dimensional stability
- Highly detailed reproduction
- Ideal viscosities
- Thixotropic
- Ideal setting time

One clear system for all your indications

Assortment:

AS40	Cavex SiliconA putty (yellow) (2 x 250 ml)
AS52	Cavex SiliconA heavy body (green) (2 x 50 ml cartridges)
AS53	Cavex SiliconA injection (blue) (2 x 50 ml cartridges)
AS70	Cavex SiliconA dispensing gun
AS94	Cavex SiliconA mixing tips injection (48 pcs.)
AS95	Cavex SiliconA mixing tips heavy body (48 pcs.)
AS96	Cavex SiliconA intra-oral tips injection Type (96 pcs.)
AS97	Cavex SiliconA intra-oral tips heavy body (96 pcs.)

Cavex Silicona: the complete, versatile range that offers the perfect solution for every challenge

Indications	Crown & bridge Sandwich	Implant Monophase	Functional Putty Wash
Putty	X		X
Heavy Body	X	X	
Injection type	X		X

YOUR IMPRESSION IS OUR SPECIALTY