

Cavex ImpreSafe



for the disinfection of imprints

A revolutionary disinfecting and cleaning fluid for alginate, polyether and silicone dental imprints. Safe, fast and effective. This agent does not contain any harmful ingredients such as aldehydes. Processing time is only 3 minutes, and this reduced disinfection time ensures that the imprint surface remains 100% intact after disinfection.

Cavex ImpreSafe

| | | |
|-------|-----------------------------|----------|
| HA053 | Cavex ImpreSafe Starter kit | 1 pce |
| HA050 | Cavex ImpreSafe Refill | 6x1000ml |

Cavex GreenClean



pH-neutral alginate and gypsum dissolvent

For the removal of alginate and gypsum residues from trays, spatulas and other instruments

Product features:

- Quick and easy
- pH-neutral, gentle on the skin
- Biodegradable, environmentally friendly
- 1 kg treats more than 500 impression trays
- Fresh smell

Cavex GreenClean

| | | |
|-------|----------------------|------|
| HA001 | Cavex GreenClean jar | 1 kg |
|-------|----------------------|------|

Cavex LC Dental Tray



light-curing tray material

Cavex LC Dental Tray is a light-curing synthetic resin for the production of individual or functional impression trays. In addition, it is suitable for other areas of application, including the production of bite templates. Cavex LC Dental Tray distinguishes itself from other brands by, amongst other things, its extra wide shape, meaning that the fuss of glueing on small pieces is a thing of the past. Furthermore, this material can be processed very easily. As standard, the tray is fit for upper jaw models; removal of the palate section, however, also makes it suitable for lower jaw models. Cavex LC Dental Tray can be cured using either UV or halogen equipment.

Cavex Light Cure Dental Tray

| | | |
|-------|------------------------------|-------|
| BP010 | Cavex Light Cure Dental Tray | 50pcs |
|-------|------------------------------|-------|

Cavex Bleach Tray Spacer



Light-curing synthetic resin

The Cavex Bleach Tray Spacer is a synthetic resin that can be easily processed. This product was specifically developed to create reservoirs in bleaching trays. It is applied directly from the handy 2ml syringe with its curved applicator tip. Ten seconds of light curing per surface is sufficient to ensure total hardening of Cavex Bleach Tray Spacer.

Cavex Bleach Tray Spacer

| | | |
|-------|--------------------------|------------------|
| VT025 | Cavex Bleach Tray Spacer | 2x2 ml + 10 tips |
|-------|--------------------------|------------------|



Cavex
Lab brochure

Cavex VacuFormer System



The Cavex VacuFormer System is a comprehensive thermal vacuum-forming system consisting of the Cavex VacuFormer, deep-draw plates and a unique gypsum composite, plus accessories.

The Cavex VacuFormer is a no-nonsense all-purpose thermal vacuum former. This rugged device features a high-performance electric vacuum motor that can cope with any conventional deep-draw application. The Cavex VacuFormer is easy to control and is protected from overheating.

The starter kit includes all the equipment and materials required to produce ready-made deep-draw applications.

Cavex VacuFormer System

| | |
|-------|---|
| VF001 | 1 Cavex VacuFormer + 18 jars of Instant Stone + deep-draw plates (50 × Splint + 25 × Bleach + 12 × Bruxism + 12 × Mouth Protector) + scissors + 6 base trays + full instructions per indication |
|-------|---|

Cavex VacuFormer System Accessories



Cavex VacuFormer System Accessories

| | | | |
|-------|---|-----------------------|-----------|
| VT020 | Cavex Mouth Protector Red/White/Blue | 3.8mm x 125mm x 125mm | 12 sheets |
| VT034 | Cavex Mouth Protector Red | 3.8mm x 125mm x 125mm | 12 sheets |
| VT035 | Cavex Mouth Protector Blue | 3.8mm x 125mm x 125mm | 12 sheets |
| VT036 | Cavex Mouth Protector Transparent | 3.8mm x 125mm x 125mm | 12 sheets |
| VT031 | Cavex Mouth Protector X-pro Green | 5mm x 125mm x 125mm | 12 sheets |
| VT032 | Cavex Mouth Protector X-pro Black | 5mm x 125mm x 125mm | 12 sheets |
| VT033 | Cavex Mouth Protector X-pro Transparent | 5mm x 125mm x 125mm | 12 sheets |
| VT021 | Cavex Splint Transparent | 0.6mm x 125mm x 125mm | 50 sheets |
| VT026 | Cavex Splint X-heavy Transparent | 1.5mm x 125mm x 125mm | 25 sheets |
| VT022 | Cavex Bleach Transparent | 1mm x 125mm x 125mm | 25 sheets |
| VT030 | Cavex Bleach Heavy Transparent | 1.5mm x 125mm x 125mm | 25 sheets |
| VT027 | Cavex Bruxism | 2mm x 125mm x 125mm | 25 sheets |

Cavex Set Up Wax & Yellow Wax



Cavex Set Up Modelling Wax

- for full and partial prostheses
- low shrinkage stress
- dimensionally stable
- plastic and elastic modelling and boxing wax

Cavex Yellow wax – soft wax

- for adjusting impression trays
- for blocking out undercuts

Cavex Set Up Was

| | | | |
|-------|----------------------|-------|------------|
| WA004 | Cavex Set Up Soft | 500g | 23 sheets |
| WA005 | Cavex Set Up Regular | 500g | 23 sheets |
| WA006 | Cavex Set Up Hard | 500g | 23 sheets |
| WA014 | Cavex Set Up Soft | 2500g | 115 sheets |
| WA015 | Cavex Set Up Regular | 2500g | 115 sheets |
| WA016 | Cavex Set Up Hard | 2500g | 115 sheets |

Cavex Yellow Wax

| | | |
|-------|------------------|-----------|
| WB005 | Cavex Yellow Wax | 20 sticks |
| WB006 | Cavex Yellow Wax | 17 sheets |

Cavex Dental Baseplates & Individual Trays



High-quality shellac plates

Cavex Baseplates

| | | |
|-------|---|--------|
| BP004 | Cavex Baseplates upper | 100pcs |
| BP054 | Cavex Baseplates lower | 100pcs |
| BP024 | Cavex Baseplates upper, double thickness; 2.70-2.80mm | 50pcs |
| BP005 | Cavex Baseplates mica-free upper | 100pcs |
| BP055 | Cavex Baseplates mica-free lower | 100pcs |

Cavex Individual Trays

| | | |
|-------|------------------------------|--------|
| BL006 | Cavex Individual Trays upper | 100pcs |
| BL056 | Cavex Individual Trays lower | 100pcs |
| BL036 | Cavex Individual upper, XL | 50pcs |
| BT009 | Handles for Individual Trays | 8pcs |
| BT019 | Handles for Individual Trays | 200pcs |



Cavex Holland BV
Fustweg 5
P.O. Box 852
2003 RW Haarlem
The Netherlands

T +31 (0)23 530 77 00
dental@cavex.nl
www.cavex.nl

Helping dental professionals and patients requires solution-based thinking.

That's why at Cavex we don't focus on separate products, but rather on high-quality systems that provide you with total solutions. Obviously, these systems comply with the highest demands in terms of quality and reliability. And it makes us proud to know that our products are used by thousands of professionals all over the world. Every day.