







Cavex Quadrant offers a carefully matched range of adhesive materials that provide the appropriate solution for every "quadrant". It includes several composites, bonding systems and a desensitizer. The Q-System is specially developed for the modern dental practice.

"Your impression is our concern"

Quadrant flow

The versatile, easy to apply flowable

-  Thixotropic - flows easy, not runny
-  Available in four different shades
-  Non sticky formula
-  Wide application range

**NOW available in shades
A2, A3, A3.5
and OA2**





CAVEX



Cavex Holland BV
P.O. Box 852
2003 RW Haarlem
The Netherlands

T +31 23 530 77 00
dental@cavex.nl
www.cavex.nl

Quadrant flow

The versatile, easy to apply flowable

Quadrant Flow is a flowable, radiopaque light curing composite. A versatile product to be used for many indications.

Quadrant Flow is easy to handle and apply. It can be placed directly on the spot, since it comes in small, convenient syringes with easy to attach delivery tips. Quadrant Flow is compatible with all modern composite systems. Undeniably, the most unique feature of Quadrant Flow, is that it is thixotropic. It flows easy but is still not runny.

Indications

Quadrant flow can be used for filling small class I cavities, fissure sealing, blocking out undercuts in inlay preparations, as an adaptive layer in deep boxes (base or liner), for filling class V cavities, for bracket and splint placement and composite repair. The OA2 shade is developed for masking purposes.

Features

- flows easy, not runny
- easy to apply
- available in four shades: A2, A3, A3.5 and OA2
- non sticky formula
- radiopaque
- compatible with other composite materials
- 3 years shelf life
- 1 ml syringes incl. delivery tips
- fluoride releasing



Assortment:

FD09	Quadrant Flow A2	1,8 g / 1 ml + 5 tips
FD10	Quadrant Flow A3	1,8 g / 1 ml + 5 tips
FD11	Quadrant Flow A3.5	1,8 g / 1 ml + 5 tips
FD12	Quadrant Flow OA2	1,8 g / 1 ml + 5 tips