

Cavex VacuFormer Laminates: Specifications

Laminate type	Splint	Bleach	Mouth Protector	Bruxism & Splint
Dimensions (inch)	5X5	5x5	5x5	5x5
Dimensions (mm)	12.7	12.7	12.7	12.7
Thickness (inch)	0.022\0.039\0.059	0.039\0.059	0.15\0.16\0.20	0.079
Thickness (mm)	0.6\1\1.5	1\1.5	3.8\4\5	2
Composition	Acrylic	Soft EVA	EVA	Acrylic