

PREMIUM RESTORATIVES

CAVEX GLASS IONOMER CEMENT

Aesthetic cement for permanent cementation
of crowns, bridges and inlays.



**YOUR
IMPRESSION
IS OUR
SPECIALTY**

CAVEX

PREMIUM RESTORATIVES

CAVEX GLASS IONOMER CEMENT

The ideal dental cement should possess several characteristics. It should be non-irritating to the dental tissue, simple to use, easy to clean up, and offer good adhesion to a variety of dental materials. Also, it should be identifiable on a radiograph (radiopaque) and provide sustained results over a prolonged period of time. Cavex developed the ideal solution to fulfill these needs:

Cavex Glass Ionomer Cement. The properties make it ideal for the cementation of crowns, bridges and inlays.

Application

Prepare the cement by gradually adding small amounts of powder to the liquid and mix to enhance absorption. Formulas of both powder and liquid are developed and optimized to attain a cement with aesthetic properties.

Features

- **Biocompatible**
- **Optimal dental tissue adhesion**
- **User-friendly processing and curing time**
- **Increased compressive and tensile strength**
- **Radiopaque**
- **No over-acidification**
- **5.5 year shelf life**



YOUR IMPRESSION IS OUR SPECIALTY

Fustweg 5 | 2031 CJ Haarlem | The Netherlands | cavex.nl

CAVEX

ZYC002EN.00