

WORKPLACE SAFETY CARD
Cavex ImpreSafe



Inside low Risk

Physical state : Liquid.

Colour : Blue.

Most important hazards :

Acute Tox. 4 H302 : Harmful if swallowed.

Skin Corr. 1B H314 : Causes severe skin burns and eye damage.

Eye Dam. 1 H318 : Causes serious eye damage.

Aquatic Acute 1 H400 : Very toxic to aquatic life.

Contains : Quaternary ammonium compounds, Benzyl-C12-18-alkyldimethyl, chlorides, isotridecanol, 3-5 EO

Hazard symbol or symbols :



Protective measures :
Not applicable.

Fire / explosion : Extinguishing media
Use extinguishing media suitable for surrounding materials.

Protective measures :
Not applicable.

INHALATION: First aid measures

INGESTION: Pure product.
Have conscious person drink several glasses of water or milk.
Aqueous, ready for use solution. Have conscious person drink several glasses of water or milk. Get medical attention / advice if you feel unwell.

Protective measures :
Avoid contact with skin. Wear suitable Nitril gloves.

SKIN CONTACT : Pure product.
After contact with skin, take off immediately all Gloves. contaminated clothing and wash immediately with plenty of soap and water. Get medical attention if irritation develops.
Aqueous, ready for use solution. Take off immediately all contaminated clothing. Wash contaminated skin with plenty of water. Get medical attention if irritation develops.

Protective measures :
Avoid contact with eyes.

EYE CONTACT :
Flush with plenty of water for at least 15 minutes, Occasionally lifting the upper and lower eyelids. Get medical attention as soon as possible and show this card to the medical personal.

Protective measures :

Wear gloves and safety glasses.



Methods of disposal :

Pure product should be regarded as hazardous waste. Aqueous, ready for use solutions should be regarded as non-hazardous after 1:3 to 1:5 dilution of the working solution. It should be disposed via the sewage system.

Storage :
Other information:

No specific storage is required.
None.

Version number : 3

Date : 22-02-2024